



Machine Learning and COVID-19

Agenda

01

The Problem



02

ML and COVID



03

Findings



04

Discuss
Opportunities




A Public Health Crisis




The virus

SARS coronavirus,
first reported Dec
2019




What we know

Developed vaccine,
effective testing




What we don't


Human responses
to virus are varied,
not fully understood



Machine Learning Solution



The Data



COVID Negative

746 CT-scans of lungs from
research papers


349 COVID positive, 397
COVID negative

Png format, image sizes
varied




COVID Positive

Building a base CNN model




1. Convolutional block
 - 2 convolutional layers
w/ 2 dropout layers
2. Convolutional block
 - 2 convolutional layers
w/ 2 dropout layers
3. Convolutional block
 - 2 convolutional layers
w/ 2 dropout layers
4. Fully connected block
 - 3 dense layers of
decreasing size

Building a base CNN model



63.1%

Applying Transfer Learning




1. VGG16 Block
2. Fully connected block
 - 3 dense layers of decreasing size
 - 2 dropout layers
3. Output layer
 - Sigmoid activation


Applying Transfer Learning

Tuning

- Implemented number of dropout layers
- Employed SGD optimizer
 - Employed high momentum and scaling learning rate
- Early stopping



Applying Transfer Learning






72.2%

Takeaways

Current model is NOT a valid diagnostic tool

ML shows promise in interpreting COVID image data

Incorporate better preprocessing, more diverse data



Thank you!

Questions?



matthewryan33@gmail.com



/maneaterrrbug




/matt-mcg-ryan

Appendix I

Layer (type)	Output Shape	Param #
conv2d_17 (Conv2D)	(None, 256, 256, 16)	448
dropout_16 (Dropout)	(None, 256, 256, 16)	0
conv2d_18 (Conv2D)	(None, 256, 256, 16)	2320
dropout_17 (Dropout)	(None, 256, 256, 16)	0
max_pooling2d_6 (MaxPooling2D)	(None, 128, 128, 16)	0
conv2d_19 (Conv2D)	(None, 128, 128, 32)	4640
dropout_18 (Dropout)	(None, 128, 128, 32)	0
conv2d_20 (Conv2D)	(None, 128, 128, 32)	9248
dropout_19 (Dropout)	(None, 128, 128, 32)	0
max_pooling2d_7 (MaxPooling2D)	(None, 64, 64, 32)	0
conv2d_21 (Conv2D)	(None, 64, 64, 64)	18496
dropout_20 (Dropout)	(None, 64, 64, 64)	0
conv2d_22 (Conv2D)	(None, 64, 64, 64)	36928
dropout_21 (Dropout)	(None, 64, 64, 64)	0
global_average_pooling2d_3 (GlobalAveragePooling2D)	(None, 64)	0
flatten_3 (Flatten)	(None, 64)	0
dense_7 (Dense)	(None, 20)	1300
dense_8 (Dense)	(None, 20)	420
dense_9 (Dense)	(None, 1)	21
=====		
Total params: 73,821		
Trainable params: 73,821		
Non-trainable params: 0		

Appendix II

Validation set accuracy by number of epochs
(base CNN)




Appendix III

Layer (type)	Output Shape	Param #
vgg16 (Functional)	(None, 8, 8, 512)	14714688
flatten (Flatten)	(None, 32768)	0
dense (Dense)	(None, 1000)	32769000
dropout (Dropout)	(None, 1000)	0
dense_1 (Dense)	(None, 500)	500500
dropout_1 (Dropout)	(None, 500)	0
dense_2 (Dense)	(None, 250)	125250
dense_3 (Dense)	(None, 1)	251
=====		
Total params: 48,109,689		
Trainable params: 33,395,001		
Non-trainable params: 14,714,688		









Appendix IV

Validation set accuracy by number of epochs
(Transfer CNN)



Appendix V

Example CT-scans



Instructions for use (free users)

In order to use this template, you must credit [Slidesgo](#) by keeping the Thanks slide.

You are allowed to:

- Modify this template.
- Use it for both personal and commercial purposes.

You are not allowed to:

- Sublicense, sell or rent any of Slidesgo Content (or a modified version of Slidesgo Content).
- Distribute this Slidesgo Template (or a modified version of this Slidesgo Template) or include it in a database or in any other product or service that offers downloadable images, icons or presentations that may be subject to distribution or resale.
- Use any of the elements that are part of this Slidesgo Template in an isolated and separated way from this Template.
- Delete the “Thanks” or “Credits” slide.
- Register any of the elements that are part of this template as a trademark or logo, or register it as a work in an intellectual property registry or similar.

For more information about editing slides, please read our FAQs or visit Slidesgo School:

<https://slidesgo.com/faqs> and <https://slidesgo.com/slidesgo-school>

Appendix